# ARAS Cranial Access Kit



IRRAS-201\_A

## **IRRAS Cranial Access Kit**



A kit designed to **better meet your cranial access needs** in a cost effective manner




## IRRAS Cranial Access Kit Key Differences

## **High Quality Components**



**Providing high-quality instruments** with every kit

### **Clear Drape**



**Clear Drape** to assist with finding mid-line

#### Larger Half Drape




**Larger Half Drape** that allows greater flexibility in creating a sterile field



## IRRAS Cranial Access Kit Competitive Analysis (1 of 4)

#### **Preparation Components**

	• •			
	Product Description	IARAS		natus.
	Double Edger Razor	<ul> <li>Image: A second s</li></ul>	✓	<ul> <li>Image: A second s</li></ul>
	DuraPrep 6ml*	<b>~</b>	✓	X
	Gauze Sponges, 4"x 4"	<b>~</b>	<ul> <li>Image: A second s</li></ul>	<ul> <li>Image: A second s</li></ul>
	Gauze Sponges, 2"x2"	X	✓	✓
	Medicine Cups, 2oz	<ul> <li>Image: A second s</li></ul>	✓	<ul> <li>Image: A second s</li></ul>
	Culture Tube with screw cap	<b>~</b>	✓	✓
	4" x 4" 6-ply Gauze	~	X	X
	Providone Iodine Solution	~	✓	~
	PVP Single Swabs	<ul> <li>Image: A second s</li></ul>	✓	X
	Povidone Iodine Ointment	X	✓	✓
	Benzoine Tincture Vial or Swab	X	X	<ul> <li>Image: A second s</li></ul>



4

## IRRAS Cranial Access Kit Competitive Analysis (2 of 4)

#### **Cranial Access Preparation**

	Product Description	IARAS		natus.
	Ruler 6"	<ul> <li>Image: A second s</li></ul>	<ul> <li>Image: A second s</li></ul>	<ul> <li>Image: A second s</li></ul>
	Marking Pen	<b>~</b>	✓	✓
	Fenestrated Drape with Barrier, 24" x 36"	X	<ul> <li>✓</li> </ul>	X
	Fenestrated Drape with Barrier, 31" x 51"	X	~	<ul> <li>✓</li> </ul>
	Fenestrated Drape with Barrier, 15" x 15"	X	~	X
	Fenestrated Drape with Barrier, 24" x 26" x 4", Clear	✓	X	X
	Film Dressing	X	X	✓
	Half Drape, 53" x 77"	<ul> <li>✓</li> </ul>	X	X
	CSR Wrap, 36" x 36"	<ul> <li>Image: A second s</li></ul>	X	X
	Absorbent Towel, 17" x 20"	<b>~</b>	X	X
	Absorbent Towel, 18" x 26"	X	✓	<b>v</b>



#### What to Remember



IRRAS' Fenestrated Drape is Clear and assists with finding mid-line easily.



IRRAS' **Half Drape** is larger and allows greater flexibility in creating a sterile field.

## **IRRAS** Cranial Access Kit **Competitive Analysis** (3 of 4)

#### **Cranial Access Components**

	Product Description	IARAS		natus.
	Lidocaine with EPI, 30ml ampule*	<ul> <li>Image: A second s</li></ul>	✓	<ul> <li>Image: A second s</li></ul>
	Sodium Chloride 10ml ampule*	<ul> <li>Image: A second s</li></ul>	X	✓
	Syringe, 10cc	X	X	<ul> <li>Image: A second s</li></ul>
	Syringe, 12cc	<ul> <li>Image: A second s</li></ul>	✓	X
	Syringe, 5cc	X	✓	X
	Needles, 18G x 1 1/2" & 25G x 5/8"	✓	<ul> <li>✓</li> </ul>	<ul> <li>✓</li> </ul>
	Needle, 23G x 1"	X	<ul> <li>Image: A second s</li></ul>	X
all a	🔨 Spinal Needle, 18G x 3 1/2"	<b>√</b>	✓	✓
	Ventricular Needle, 12G x 5 1/2"	X	<ul> <li>Image: A second s</li></ul>	X

	Product Description	IARAS		natus
	Scalpel, #11 & #15	<ul> <li>Image: A second s</li></ul>	<ul> <li>✓</li> </ul>	<ul> <li>Image: A second s</li></ul>
	Retractor, Blunt	<ul> <li>✓</li> </ul>	✓	✓
Total .	Hand Crank Drill	<ul> <li>Image: A second s</li></ul>	✓	✓
M	Aosquito Forceps, curved	<ul> <li>✓</li> </ul>	✓	<ul> <li>✓</li> </ul>

What to Remember



IRRAS provides highquality instruments with every kit.



IRRAS Kit - IC40-D, Integra Kit - INSHITHRZN (https://bit.ly/2xsrtEN), Natus Kit - HITHRZR (https://bit.ly/33KpoQV) \* Specific components included in the IRRAS IC40-D kit

## IRRAS Cranial Access Kit Competitive Analysis (4 of 4)

#### **Drill Bits**

Product Description	IARAS		natus.
2.7mm Drill Bit*	<ul> <li>Image: A second s</li></ul>	X	<b>v</b>
4mm Drill Bit	✓	1	X
5.3mm Drill Bit	X	<ul> <li>Image: A second s</li></ul>	<ul> <li>Image: A second s</li></ul>
5.8mm Drill Bit*	✓	X	X
6.35mm Drill Bit	X	<ul> <li>Image: A second s</li></ul>	X
7.1mm Drill Bit	~	X	X

## Wound Closing/Dressing

Product Description	IARAS		natus
Tegaderm	<b>v</b>	X	X
2.0 Silk Sutures**	<b>~</b>	X	✓
3.0 Nylon Suture**	<ul> <li>Image: A second s</li></ul>	<ul> <li>Image: A second s</li></ul>	✓
Needle Holder, Serrated	<ul> <li>✓</li> </ul>	<	<ul> <li>✓</li> </ul>
Adson Forceps	<ul> <li>Image: A second s</li></ul>	✓	<ul> <li>Image: A second s</li></ul>
Adson Forceps with Teeth	<ul> <li>✓</li> </ul>	✓	✓
Suture Scissors	<b>v</b>	<ul> <li>✓</li> </ul>	<ul> <li>Image: A second s</li></ul>



#### What to Remember



IRRAS offers a variety of different sized drill bits to customize each treatment.

IRRAS Kit - IC40-D, Integra Kit – INSHITHRZN (<u>https://bit.ly/2xsrtEN</u>) , Natus Kit – HITHRZR (<u>https://bit.ly/33KpoQV</u>)

\*Components included in the IRRAS IC42 kit \*\* Specific components included in the IRRAS IC40-D kit



## **FRAS** is bringing needed innovation to neurocritical care by providing a cranial access kit designed to meet your needs

## A simplified solution for cranial access

## **High Quality Components**



**Providing high-quality instruments** with every kit





**Clear Drape** to assist with finding mid-line

### Larger Half Drape



Larger Half Drape that allows greater flexibility in creating a sterile field



## **FRAS** Needed Innovation for Neurocritical Care Offering a full product portfolio from the IRRAflow and Hummingbird product lines



Hummingbird

## Product family portfolio

Product family portfolio



**IRRA***flow* Intelligent **Digital Pump** 

> Enables interaction between control unit and catheter



Controlled management of ICP and drainage via touchscreen console

**Cranial Access Kit** 

Essential tools for cranial access, used in almost every neurosurgical operation

**ICP** Monitoring **Hummingbird Solo** 

**ICP** Monitoring . . . . . . . . . . . . . . . .

**Combines brain** access and parenchymal ICP measurement

Parenchyma = Brian tissue

**Hummingbird Quad ICP** Monitoring Combines brain access. ICP measurement, drainage, and two additional lumens

## **Hummingbird ICP Control Module**

Connects Hummingbird disposables to ~10 different **Hospital Monitors**